

Jazzy
Carbon

Ultra-Lite



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

Weighing only 14,5 kg and crafted from durable carbon fibre, the Pride Jazzy Carbon Ultra Lite offers a sleek, modern design that blends style with functionality.



One lever fold function



Fold away footplate



Stretch mesh backrest for light-weight extra comfort

Specifications

Weight capacity	125 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 20 km
Turning radius	900 mm
Length	1000 mm
Width	620 mm
Drive wheels	12" pneumatic or solid
Caster wheels	8" solid
Weight without batteries	14,5 kg
Batteries	Type: Lithium ion Size: 10 Ah

Features

- At just 14,5 kg, the Pride Jazzy Carbon Ultra Lite is easy to transport, fold, and store.
- Travel at speeds of up to 6,4 km/h with ease, whether indoors or on flat outdoor pavements.
- With an onboard USB charger, you can keep your devices powered up while on the move.
- The one-touch folding system makes it quick and easy to fold the chair, perfect for travel and storage.
- Equipped with lithium batteries that are safe for air travel, making your next adventure hassle-free.

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



PRIDE
MOBILITY
LIVE YOUR BEST®

Pride Mobility Products Europe
De Zwaan 3
1601 MS Enkhuizen
The Netherlands
Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
pridemobility.eu · info@prideeurope.nl



EU_SSP0057/Rev0/27NOV24